**ZAR**  
**ACCOM**  
L E N S E S

CE 1984

TR

EN

**ZAR**  
**ACCOM**  
L E N S E S

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ÜRETİCİ

ANADOLU TIP TEKNOLOJİLERİ

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**ZAR**  
**ACCOM**  
L E N S E S

**AS-05/2022-RV 02**

**Foldable Aspheric Intraocular Lens Instruction For Use**

Zaraccomm Aspheric (AS60125),

Zaraccomm Aspheric L (AS60130), Zaraccomm Aspheric (AS60130L)

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
**DESCRIPTION:** Zaraccom foldable aspheric intraocular lenses (IOLs) are designed for cataract surgery as an anterior surface aspheric optic device to replace the crystalline lens and to reduce the spherical aberration in the human eye. Foldable IOLs consist of hydrophobic acrylic and are manufactured with bonded UV-absorber, are sterile, 1-piece and suitable for posterior chamber implantation. They have a biconvex optic with supporting haptics.

**Table1. Physical Characteristic of IOLs**

Physical Characteristic	Description
Models	AS60125, AS60130, AS60130L
Material	Hydrophobic Acrylic (Acrylate-Methacrylate Copolymer)
Index of Refraction	1.51
Optical Design	Aspheric, Bi Convex, square edge
YAG Laser	Stable
Sterilization Type	Ethylene Oxide
Power	0 Dpt through +40 Dpt (0.50 Dpt increments)
Estimated A. Constant (Calculated theoretically)	118.4
Haptic Angle	0°
Haptic Configuration	Modified L-Loop
Water Content	<1%

**Figure 1: Design of Zaraccom IOL Models**

### BIOCOMPATIBILITY TESTING

Potential patient safety risks to the material(s) of this device were evaluated through nonclinical physicochemical characterization and biocompatibility testing in accordance with international standards applicable to IOL devices. Nonclinical testing demonstrated no safety concerns for local or systemic toxicity, that the IOL material was physically and optically stable, and that there were no leachable substances arising from the manufacturing process (including sterilization) or device material(s) that posed a safety risk. The device possesses an acceptable patient safety profile when used in accordance with the Instruction for Use (IFU) for its intended clinical purpose as an ocular implant device.

## INTENDED USE

Zaraccom posterior chamber intraocular lenses are intended for use by a trained ophthalmic surgeon. The IOL is intended to be positioned in the posterior chamber of the eye, replacing the natural crystalline lens. This position allows the lens to function as a refractive medium in the correction of aphakia. These lenses are intended for placement in the capsular bag.

## INDICATIONS

Zaraccom foldable posterior chamber IOLs are indicated for primary implantation for the visual correction of aphakia in adult patients in whom a cataractous lens has been removed.

## CONTRAINDICATIONS

To date there are no absolute contraindications to Intraocular Lenses implantation. There are no known contraindications with the use of Zaraccom IOL when used as recommended.

## RESIDUAL RISKS AND UNDESIRABLE EFFECTS

As with any surgery, there is risk involved whether or not an IOL is implanted. Potential complications accompanying cataract or implant surgery are not limited to the abovementioned complications.

Cataract surgery, with or without lens implantation might be associated with:

- Ocular inflammation
- Hemorrhage
- Intraocular pressure elevation
- Post-operative infection
- Retinal breaks and detachment
- Cystoid macular edema
- Corneal edema
- Posterior capsule opacity

Complications related with Intraocular Lens implantation:

- Capsular rupture
- Vitreous loss
- Lens decentration and luxation
- Wrong calculation of IOL power
- Damage of IOL during implantation

## WARNINGS AND PRECAUTIONS

- A high level of surgical skill is required for intraocular lens implantation. The surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on intraocular lens implantation before attempting to implant intraocular lenses.
- Prior to surgery, prospective patients should be informed of the possible risks and benefits associated with this IOL as well as the risks and benefits associated with cataract surgery.
- Do not use the intraocular lens in case the sterile package was opened or damaged or in any case of doubt.
- Do not use the intraocular lens after expiration date.
- Do not resterilize this intraocular lens by any method.
- Do not reuse the lens. Single use only. The lens must be used once only for a single patient. Reuse of this single-use device may result in serious injury, such as endophthalmitis.
- Do not store at under 5°C (41°F) and above 35°C (95°F). The manufacturer recommends the storage and usage of lenses at room temperature (This allows the intraocular lens to open more slowly in the capsule). Do not expose the lens to sunlight.
- Do not soak the lens in solutions other than balanced saline solutions or equivalents.
- This lens is used safely with viscoelastic materials (OVD).

- Handle the lens carefully to avoid damage to lens surfaces or haptics.
- Do not attempt to reshape haptics in any way.
- After implantation the lens position must be Reversed-S.
- The long-term effects of intraocular lens implantation have not been determined. Therefore, physicians should continue to monitor patients postoperatively on a regular basis.
- Zaraccom intraocular lenses are posterior chamber lenses. The safety and effectiveness of a posterior chamber lens, if placed in the anterior chamber, has not been established.
- It is recommended that OVD be removed from the eye at the close of surgery with emphasis on the space between the posterior capsule and lens. This may be accomplished by gently depressing the IOL optic posteriorly with the irrigation/aspiration (I/A) tip while using standard I/A techniques to remove the OVD from the eye. This should force any trapped OVD anteriorly where it can be easily aspirated.

## SUGGESTED A-CONSTANT:

The suggested A-constant indicated on the outer label is presented as a guideline and is a starting point for implant power calculations. It is recommended that you develop your own constant appropriate for you based on clinical experience with the particular lens models, surgical techniques, measuring equipment, and postoperative results. If additional information on lens power calculation is needed, please contact with the manufacturer. The contact data are given at manufacturer details.

## CALCULATION OF LENS POWER:

Accurate biometry is essential for successful visual outcomes. Preoperative calculation of required lens power for these posterior chamber intraocular lenses should be determined by the surgeon's experience, preference, and intended lens placement. A reference A-Constant value is located on the outer label. Lens power calculation methods are described in the following references:

Hoffer, K.J. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. 19:700-712, 1993.

Holladay, J.T., et al. A Three-part system for Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. 14:17-24, 1988.

Holladay, J.T., et al. Standardizing Constants for Ultrasonic Biometry, Keratometry, and IOL Power Calculations. J. Cataract Refract. Surg. 23:1356-1370, 1997.

Retzlaff, J.A., Sanders, D.R., and Kraff, M. Lens Implant Power Calculation, 3rd ed., Slack, Inc., Thorofare, N.J., 1990.

Haigis W: The Haigis Formula, In: Intraocular lens power calculations. H. John Sammas (eds) Slack Incorporated, Thorofare NJ, USA, pp 39-57, 2004

## INSTRUCTIONS FOR USE:

1. Prior to implanting, examine the lens package for proper lens model, dioptic power, and expiration date.
2. Open the package and remove the lens in a sterile environment. The lens is provided in a sterile plastic container. To remove the lens, internal box cover is slowly pulled backward as shown by the arrow direction. The protective cover lifts and the lens is taken out of the internal box with forceps. The lens is placed in the proper position for the implantation. The lens must be taken and implanted carefully, in order not to scratch the optic surface by an instrument.
3. Rinse the lens thoroughly using sterile balanced salt solution.
4. To minimize the occurrence of marks on the lens due to folding, all instrumentation has to be scrupulously clean. The manufacturer recommends using forceps with round edges and smooth surfaces.
5. Handle the IOL by the haptics only when removing the lens from the case. DO NOT grasp the optical area with forceps.
6. The lenses are carefully checked and inspected by manufacturer to assure a high quality product. The lenses are to be studied carefully under surgical operation microscope before implantation. If a defect or deformation is noted or suspected, the lens should be returned to the manufacturer.
7. Implantation has to be performed on aphakic adult patients by a trained ophthalmic surgeon.

8. Implant the IOL with the most appropriate surgical procedure for the patient. Surgeons should verify that appropriate instrumentation is available prior to surgery. Refer to ANY OTHER DEVICES AND PRODUCTS WHICH ARE INTENDED TO BE USED IN COMBINATION WITH THE DEVICE for qualified or compatible products.

### **ANY OTHER DEVICES AND PRODUCTS WHICH ARE INTENDED TO BE USED IN COMBINATION WITH THE DEVICE**

The foldable lenses can be folded by various folding forceps designed for foldable lenses and it can be implanted also by implantation forceps and cartridge-injector systems. Posterior chamber lenses are designed for implantation in the capsular bag by injector method with a minimum incision size of 2.2 mm. It is recommended to use cartridge systems that validated and approved by manufacturer. The use of an unqualified combination may cause damage to the lens and potential complications during the implantation process (e.g., IOL decentration, tilt, or dislocation, spatial distortions, visual disturbances, or decreased/blurred vision). In addition, if the instructions for use are not followed and/or used by persons other than a trained ophthalmic surgeon, "Abnormal Use" will occur, and the following hazardous situations may be encountered. Although this situation has no effect on the patient, it causes the intraocular lens not to be used in the operation.

Abnormal Use Cases;

- The intraocular lens gets stuck in the cartridge injector system,
- Haptic rupture, breakage,
- Scratches, cracks on the optical surface,
- Uncontrolled protrusion of the intraocular lens from the injector.

Qualified cartridge injector systems that can be used with Zaraccom foldable IOLs are listed at IFU ANNEX document.

### **MAGNETIC RESONANCE COMPATIBILITY**

Zaraccom foldable hydrophobic IOL is magnetic resonance (MR) Safe. The IOL consists of acrylate/methacrylate copolymer material, which is a non-conducting, non-metallic, non-magnetic material that poses no known hazards in all magnetic resonance imaging environments.

### **LIFETIME OF THE IOL**

Based upon material characterization of the Zaraccom IOL material, the IOL is expected to be stable indefinitely over the lifetime of the patient.

## PACKAGING AND EXPIRATION DATE:

Zaraccom lenses are supplied in dry, in a double sterile package terminally sterilized with ethylene oxide, and must be opened only under aseptic conditions. The packaging sterility is guaranteed until expiry date except in case of damaging or opening of the package. The use-by date is clearly indicated on the outer box label of the carton. Any intraocular lens held after the use-by date should be returned to Anadolu Tip Teknolojileri A.Ş.

## LIABILITY:

The manufacturer's liability covers the design and production of these intraocular lenses; it shall be incurred in no way in case of accidents resulting of the use of these lenses.

## SERIOUS INCIDENT REPORTING:

Any serious incident that may reasonably be regarded as device related should be reported to;

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By Phone: +90 346 218 14 18


By Email: [mail@anadolutip.com.tr](mailto:mail@anadolutip.com.tr)

Website: <http://www.anadolutip.com.tr>

Each IOL is identified by a serial number which provides traceability, and this information should be given to Anadolu Tip Teknolojileri A.Ş.

**NOTE:** In Europe, these serious incidents must also be reported to the competent authority for medical devices of the appropriate State.

## SYMBOLS USED ON LABELING:

Symbol	Explanation	Symbol	Explanation
	Sterilized using ethylene oxide.		Use by date.
	Temperature limitation: lower temperature limit is 5°C (41°F) and upper temperature limit is 35°C (95°F)		The product conforms to European Medical Directive 93/42/EEC and meets applicable health, safety and environmental requirements. If the mark is accompanied by a number, conformity is verified by the indicated notified body.
	Do not re-sterilize.		Do not expose to sunlight.
	Re-use is not allowed		Not to be used in case package is damaged
	Manufacturer		Read usage instructions

## IFU ANNEX Document

### Suitable cartridge and injector systems are as below;

Product : Top Loaded Cartridge Injector System

Manufacturer : ODC Inc. / FRANCE

Model : Ergotouch FLY 2.2 ve Ergotouch FLY 2.4

Recommended power range for Ergotouch FLY 2.2 :

Up to 25,00 Dpt Zaraccom IOLs

Recommended power range for Ergotouch FLY 2.4 :

Over 25,00 Dpt Zaraccom IOLs

The instructions for use are summarized, please refer to the device's own IFU for detailed information.

### Top Loaded Cartridge Injector System Instruction for Use

1. Open the blister pack using standard sterile procedures. Place contents onto a sterile field.
2. Fill the cartridge tunnel of the Top Loaded Cartridge Injector System and the loading chamber of the cartridge sufficiently with validated viscoelastic solution. Balanced salt solution should not be used as the sole lubricant.
3. Place the lens symmetrically on the central hinge of the open

cartridge. Asymmetric insertion of the lens can lead to rotation of the lens in the injector.

4. Press the lens into both guide rails with rounded sterile tweezers and gently move the lens back and forth horizontally in order to ensure that it can move freely.

5. Ensure that the lens is correctly aligned and make sure that neither the edge of the lens nor the haptic are wedged in.

6. Press the sides of the cartridge together until the click-lock mechanism engages. The audible click demonstrates the correct closure of the cartridge. A final visual check is recommended after closure to ensure haptic is not trapped.

7. Insert the cartridge into the front end of the injector. Push the cartridge to the very front position.

8. Carefully push the plunger forward and ensure that the silicone tip correctly enters the loading chamber. If the silicone tip cannot be introduced into the loading chamber, retract the plunger to the starting position and align the silicone tip straight, using sterile tweezers. Continue to push the cushion until the inner spring begins to press together. Pull the plunger back a few millimeters and then push it forward again. This step ensures that the lens is correctly grasped.

9. In order to avoid excess viscoelastic entering the eye, gently push the plunger until the lens is at the front of the cartridge, releasing the excess viscoelastic well away from the eye before inserting the injector into the incision.

10. Guide the tip of the cartridge through the incision and push it across the iris to the near edge of the pupil.
11. Slowly press the plunger forward to support the correct delivery of the lens into the eye and simultaneously withdraw the instrument from the eye. To avoid swelling of the silicone tip during ejection from the cartridge, only push the plunger until the lens has emerged completely.
12. Position the lens carefully and if necessary rotate it with the aid of the appropriate positioning hook.
13. Remove the viscoelastic material from the eye and the lens, with standard irrigation and aspiration techniques.

### **Product: Back Loaded Cartridge Injector System**

**Manufacturer** : AST Products Inc. / USA

**Model** : BIOLI-D, BIOLI-C

**Recommended power range for BIOLI-C:**

Up to 25.00 Dpt Zaraccom IOLs

**Recommended power range for BIOLI-D:**

Over 25.00 Dpt Zaraccom IOLs

The instructions for use are summarized, please refer to the device's own IFU for detailed information.

### **Back Loaded Cartridge Injector System Instruction for Use**

1. Inspect the blister packaging for any damage. If damage is observed or suspected, use another injector set. Open the blister packaging using standard sterile procedures and transfer the injector and cartridge to the sterile field.
2. Apply viscoelastic solution to the Back Loaded Cartridge Injector System. Ensure that viscoelastic solution is at operating room temperature before filling the cartridge, immediately prior to loading and delivery of the lens.
3. Using holding forceps, grasp the lens by the optic edge and gently place the lens anterior side up into the rear of the cartridge. The lens should be inserted until the optic is a little more than half-way inside the cartridge. Use the holding forceps to gently push down on the lens. Verifying that the lens is on the bottom surface of the cartridge.
4. Using holding forceps, grasp the trailing haptic, and gently fold the haptic onto the anterior side of the optic.
5. Slowly push the optic edge to position the lens as far into the cartridge as the forceps will permit, while ensuring the lens remains on the bottom surface of the cartridge and the trailing haptic remains on top of the optic.
6. Be careful not to scratch the IOL during this process. Remove the forceps after pushing the lens, and immediately.
7. Insert the cartridge into the handpiece by sliding the cartridge backward into the handpiece slot.

8. Advance the plunger in one slow motion. The plunger should make initial contact with the cartridge at the ramp.
9. Slowly push the plunger forward to advance and fold the optic. Verify that the lens moves forward at the same rate as the plunger.
10. Continue to push the plunger slowly until the spring engages. Slightly pull back the plunger about 1 mm before continuing to push the plunger until the lens has been safely delivered.
11. As the lens moves through the cartridge, verify that the plunger tip is behind the optic and the leading haptic is looped in front of the optic. The leading haptic should become looped back between the folded halves of the optic. If the leading haptic is not looped back, retract the plunger and reload the lens again.
12. Insert the tip of the cartridge through the incision. Position the tip at the anterior capsule opening. Advance the lens only until the optic exits the nozzle.
13. Withdraw the cartridge nozzle once the lens has exited the cartridge nozzle.
14. Place and center the lens using a suitable positioning instrument.
15. Discard the injector and cartridge accordingly.

## Precautions

- The Top/Back Loaded Cartridge Injector System should be used by trained and qualified professionals only. The result of surgical procedure depends on skill and abilities used with the device and other products during the surgery
- This product is supplied sterile and should not be used if inner packaging is soaked, opened or damaged
- All contents are for single use and should never be re-sterilized
- Do not use the cartridge if tip has been broken or damaged
- Always use viscoelastic when loading the foldable intraocular lens into the cartridge. The lens should be injected immediately after being introduced into the cartridge as viscoelastic solution may lose its lubricating properties when exposed to air for a long period of time
- Please follow up carefully the loading instructions
- In case the device falls on the floor or touch any other unsterile product, discard it according to standard rules for residues in health establishments

If the above precautions are not followed, foreseeable risks or complications could cause:

- Improper handling, loading and use of the injector can cause damage to the lens and patient such as tissue damage, inflammation, contamination
- Contamination during use could generate infection / contamination of the patient
- Re-use of the Single-use Injector and/or cartridge could generate contamination between patients
- There are risks such as the intraocular lens getting stuck in the cartridge injector system, haptic breakage, scratches and cracks on the optical surface, uncontrolled removal of the intraocular lens from the injector.